

# Duo 200

*The Duo 200, for electrotherapy.*

The Duo 200 : a 'must' for your practice. Portable and supremely easy to use, with a large backlight screen and two independent channels that allow two different types of electrotherapy currents to be applied at the same time. Pushbuttons for the most important functions. 21 types of electrotherapy current, and simultaneous treatments. Han Stim, the new pain control therapy, and treatment for an overactive bladder have been added to the therapeutic objectives.



Article number: 300 398

## Characteristics

- Therapies :
  - Electro therapy, 2 and 4 poles, 2 channels completely independent
  - Simultaneous therapy (2 different indications treated simultaneously by using electro)
- Currents :
  - Current forms: see on the next page
- Functionalities :
  - Objectives: 19
  - Indication: 85
  - Diagnostic: 6
  - Free memory: 20

Can be connected to the Vaco 500

## Technical specifications

Languages	: English, German, French, Dutch, Spanish, Italian, Portuguese and Polish
Mains voltages	: 100- 240-VAC, 50/60 Hz 50/60 HZ +/- 10 %
Maximum power -in operation	: 85 VA
Dimensions (b x h x d)	: 266 x 100 x 275 mm
Weight	: ca 3,650 kg
Safety class	: Class I
Isolationtype	: Type BF
Classification according to MDD	: Class IIa
Safety examination	: CE from conformity with Directives MDD 93/42/EEC

## Standard accessories

100 689	Mains lead
108 725	Electrode cable, 2-pole (2x)
108 919	Test plug, F/F, 4 mm
109 959	Electrode, rubber, 6 x 8 cm, per 2 (2x)
100 658	Chamex bag, 6 x 8 cm, per 4
108 935	Fixing strap, elastic, 5 x 60 cm (4x)
115 684	Visual Analogue Scale (VAS-Score)
117 276	Sets of photos Duo 200 electrodes placement

## Manuals

323 011	Safety Instructions
117 268	Treatment Guide Dutch
117 269	Treatment Guide French
117 271	Treatment Guide English
117 270	Treatment Guide German
327 382	Treatment Guide Spanish

## Optional accessories






302 955	Carrying case 200 series
114 142	Pen electrode with sponge, Ø 15 mm
109 944	Sponge for pen electrode, per 10
329 978	Vaginal probe Novatys Gold
330 594	Vaginal probe V2B+
330 572	Vaginal probe Optima 3
330 583	Vaginal probe Perisize 4+
329 989	Anal probe Analia
330 561	Anal probe Analys+
112 166	Stimulation probe, rectal
326 799	Electrode, adhesive, Ø 3 cm, per 4
326 810	Electrode, adhesive, 2.5 x 5 cm, per 4
326 821	Electrode, adhesive, 5 x 5 cm, per 4
326 832	Electrode, adhesive, 5 x 10 cm, per 4
328 504	Adhesive electrode kit
113 334	Adapter for adhesive electrode 4 > 2 mm
109 958	Electrode, rubber, 4 x 6 cm, per 2
109 960	Electrode, rubber, 8 x 12 cm, per 2
108 934	Fixing strap, elastic, 5 x 30 cm
108 936	Fixing strap, elastic, 5 x 120 cm
100 657	Chamex bag, 4 x 6 cm, per 4
100 659	Chamex bag, 8 x 12 cm, per 4
117 210	Mobil 700, trolley, 2 shelves, 1 drawer
117 220	Extra shelf + drawer for Mobil 700
117 215	Bottle support for Mobil 700








# Duo 200 : Current forms






Unidirectional Currents		
Medium frequency rectangular	Rectangular pulse	Triangular pulse
2-5 current (Ultra Reiz)		

Diadynamic Currents		
MF	DF	CP
CP-ISO	LP	LP-ISO

Tens Currents		
		
Conventional Tens example: Symmetrical	Low Frequency Tens example: Symmetrical	Random frequency Tens example: symmetrical
		
Burst Tens example: symmetrical	Han Stim	

NMES Currents		
		
Rectangular Surge	Triangular Surge	Biphasic Surge example: symmetrical
		
Biphasic Surge example: a-symmetrical	Intrapulse Interval Surge Symmetrical	2-pole Medium Frequency Surge
		
4-pole: Interferential surge		

### Interferential Currents

		
2-pole Medium frequency	4-pole Interferential	4-pole Interferential with Rotating Vector